



Neonatal · Pediatric · Adult

BIYOVENT


Intensive Care Type
Mechanical Ventilator Device



biosysmed.com



Breathe the Future ●●●



biosysmed.com



Our Work, The Source of Life

Breathing for health, health for technology

As Biosys, since 2012, we have been providing solutions to gaps in the healthcare sector with experienced engineers and expert doctors with extensive R&D experience, starting our journey in Ankara. In our working process, we advance with the awareness that access to healthcare services is a fundamental social right and we develop our technologies with the consciousness of delivering this right to every individual.

By developing next-generation medical technologies, we turn our ideas into reality at our R&D center to meet both national and global needs in the healthcare field. We know that even the smallest task we undertake can impact an individual's health somewhere in the world, and for this reason, we carry out meticulous and in-depth work to ensure that every detail of our products is designed flawlessly. Thanks to our effective and safe production approach, we produce next-generation medical technologies with high quality and reliability, becoming a beacon of hope for healthy living.

In the technologies we develop to meet the needs and demands of the healthcare sector, we do not rely solely on our engineering expertise. Instead, we engage in a multidisciplinary approach shaped by the input of doctors at every step we take. Thus, we create medical devices that not only have high technological threshold values but also fully meet the needs, encompassing the best in design, software, and technology. With this distinctive approach, we continue to move forward confidently, aiming to touch the health of every individual and become a gateway to healthy living worldwide.




18000+
Devices



40+
Countries



20
Delivery in less
than twenty days





DEFINITION OF **PERSONALIZED TREATMENT**

Optimal Treatment for Every Patient's Unique Needs

Each patient's condition and needs are different. With this awareness, Biyovent offers a customizable structure based on the patient's needs. This approach provides advantages for both the patient and the doctor, ensuring the most effective treatment is applied.

Biyovent, designed with every detail considered for health, helps to apply the most suitable treatment method for pediatric, adult, and neonatal patients. By recording values related to the patient's condition, changing clinical needs can be monitored. This ensures that the most appropriate and accurate treatment is provided to each patient.

From a usability standpoint, Biyovent provides healthcare providers with access to various necessary information while enabling the application of the required treatment method for patients. It offers the flexibility to choose between Non-Invasive, Invasive ventilation, N-CPAP, and HFO therapy types. Biyovent establishes a secure infrastructure for treatment with advanced software features, allowing adjustments such as circuit length, patient-specific circuit thickness selection, leak trigger compensation, and adjustable alarm limits. It elevates safety measures to the highest standards.





INNOVATIVE DESIGN, **COMPREHENSIVE SOLUTIONS**

Design Suitable for Quick and Easy Setup

Biyovent, with its three-part structure, can be set up quickly and easily. Its ergonomic design, consisting of base, support, and main body parts connected to the screen, facilitates ease of installation.

Biyovent simplifies solution-oriented work with a streamlined and quick installation system that goes beyond complex configurations.





Invasive and Non-Invasive Ventilation



Integrated Nebulizer



High Flow Oxygen Therapy and N-CPAP
(Nasal Continuous Positive Airway Pressure)



Suitable for Pediatric, Adult, and
Neonatal Patients



Smart Ventilation Modes



Leak and Trigger Compensation



- Multi-use autoclavable expiratory valve
- Nebulizer output



- 15 inch high-resolution color display
- Sensitive touch
- Easy transition between modes
- Easily adjustable mode parameters
- Capability to monitor more than 30 parameters simultaneously

- Wheels with increased shock absorption
- Locking mechanism
- Safe design
- Ease of use




- Power supply for humidifier
- External nebulizer input
- Emergency valve






ADVANCED SOFTWARE, EASY USE



With an easy-to-use interface design, patient information can be entered effortlessly and allows customization of parameters tailored to the patient's needs.



It allows easy calibration in 3 steps through on-screen guidance.

USER CENTERED, SCREEN SYSTEM


Biyovent, advanced software features and a user-centered interface design, simultaneously displays 3 graphs, 3 loops, and lung dynamic status. Its visually supported interface design facilitates easy access to information, enhancing learning and application. It allows access to more than 30 parameters on the same screen, and its 15-inch high-resolution display presents information clearly and in detail.



It allows seeing the dynamic status of the lungs with a lung diagram that offers smart respiratory monitoring feature.



Monitoring parameters through a single screen allows simultaneous monitoring of more than 30 parameters.





ENHANCE YOUR **BIYOVENT EXPERIENCE**

BioAqua humidifier aligns with invasive and non-invasive mechanical ventilation. It can also be used during a HFO therapy provided by a ventilator. BioAqua has invasive, non-invasive and HFO treatment modes provided with different temperatures accordingly.

The special exterior design of the design enables monitorization of parameters easily when the BioAqua is attached to the ventilator.

The device provides the nebulizer feature with standard configuration of humidification. With different modes of nebulizer flow rate, accurate treatment is achieved.



Ventilation Modes Suitable for Every Age and Condition

Pressure Targeted Modes	Volume Targeted Modes	Spontaneous Modes
P-ACV	V-ACV	SPN-PS
P-SIMV+PS	V-ACV(PRVC)	SPN-VS
P-CMV	V-CMV	Smart Modes
P-PSV	V-SIMV+PS	SAV
P-Bilevel	V-SIMV(PRVC)+P	
APRV		



TECHNICAL SPECIFICATIONS

Monitoring Parameters			
Ppeak (cmH2O)	Compliance (ml/cmH2O)	Peep (cmH2O)	RC Constant (sec)
Minute Volume (Lt/Min.)	Tidal Volume (ml)	Port (cmH2O)	Leakage(%)
Respiratory Rate (RR/Min.)	Pplato (cmH2O)t	Ti (sec)	Te(sec)
I:E Rate	Oxygen (%)	WOB (joule)	WOB/Lt (joule/Lt)
Inspiration Peak Pressure (Lt/Min)	Expiration Peak Pressure (Lt/Min)	Vexpiraton (mL)	Vinspiration (mL)
Vresidual (ml)	Exhalation S(Spon) (SS/min)	Leakage Constant(Ltr/m)/cmH2O)	Minute Volume (Spn) (Lt/min)
Elastance (cmH2O/Ltr)	OtoPEEP (cmH2O)	P0.1 (cmH2O/ 0.1s)	P_nif (cmH2O)
RSBI (rr/Min)Lt	PTP (cmH2OXsec)	rcExp (sec)	Expiration Volume/min (Lt/min)
Leakage Constant (Lt/ml/cmH2O)	Instantaneous Tidal Volume (mL)		



Patient Compliance and Performance		Pressure Source Features
Neonatal	Adult & Pediatric	O2 Pressure: 2.5 - 7 Bar
Respiratory rate: 1-150/min	Respiratory rate: 1-100/min	Air Pressure: 2.5 - 7 Bar
T inspiration: 0.1-10 second	T inspiration: 0.1-10 second	Automatic changeover and alarm at the end of supply
Tidal volume: 2 -300ml	Tidal volume: 100-3000ml	Medical compressor/Ability to operate with a regulator
Flow: 1-60 lt/min	Flow: 2-120 lt/min.	Electrical Specifications
		Batttery Life: 2 hours + 8 hours (Optional)
		Mains Voltage: 180 - 264 VAC
		Power Consumption: 47 - 63 Hz 100W
Working Features		Digital Interfaces
Inspiration pressure: 2-100 mbar	Inspiration Time: 0.1-10 second	USB: 4 pieces / COM: 2 pieces
Peep pressure: 1-50 cmH2O	Breathing Rate: (c) 1-150/min (a) 1-100/min	HDMI: 1 piece / LAN: 2 pieces
Flow: (c) 1-60 lt/min (a) 2-120 lt/min (c): Child, (a): Adult	O2 Mixture: 21-100%	Display Port: 1 piece
Spontane Pressure Support: 0-100 mBar	I/E Rate: 1:10 -10:1	Mic. In: 1 piece
		Dimensions
		B: 150 cm / D: 58 cm / G: 58 cm
		Weight: 55 Kg
		Screen Movement: Right or Left 150° Up or Down 15°
		Screen: Fully touchscreen



BIYOVENT

Intensive Care Type Mechanical Ventilator Device



Üniversiteler Mahallesi İhsan Doğramacı Bulvarı No: 23/C
ODTÜ Teknokent 06800 Çankaya/ANKARA

+90 850 800 62 97

biosysmed.com

